

NIDEK *Medical*

THE MAX 30



Model Part Number **3005**
230 Volt / 50 Hz

Model Part Number **3010**
230 Volt / 60 Hz

Specifications

Average Power Consumption **12 Amps**

Oxygen Concentration **93% (+ / -3%)
at 30 LPM**

Liter Flow **2 to 30 LPM**

Outlet Pressure **50 PSIG (3.4 BAR)**

Weight **200 lbs / 90 kg**

Dimensions **24"W x 20"D x 45.5"H
61 cm x 51 cm x 116 cm**

Flow Meters and External Storage Tanks
available upon request --
Please contact your sales rep!

Applications

Cylinder Filling Stations

Anesthesia Machines

Small Hospitals

Laboratories

Veterinary Clinics

Recovery Ventilators

Oxygen Therapy


CE
0413

2010-9806 REV B

THE MAX 30

Configurations

Simple Hospital Piping System



Small Clinic System

