



NidekMedical
PRODUCTS

.....Mark 5
Nuvo
Lite & Lite 3

CE
2862

Specifications

.....Mark 5
Nuvo
Lite



Model 525	Model 925
115 V / 60 Hz	230 V / 50 Hz
330 Watts	300 Watts

.....Mark 5
Nuvo
Lite3



Model 535	Model 935
115 V / 60 Hz	230 V / 50 Hz
210 Watts	180 Watts



90% (+5.5% / -3%)
at 3 and 5 LPM



0.125 to 3 LPM
0.125 to 5 LPM



< 48 dBA



7.0 PSI



32 lbs
14.5 kg



14" W x 9" D x 23" H
36 cm x 23 cm x 58.5 cm

2010-8422 Rev C

nidekmedical.com